





DUS 8

Digital Ultrasonic Diagnostic Imaging System

DUS 8, EDAN's new member of the all digital ultrasonic diagnostic imaging systems, not only endeavors to offer the most competitive price, but also focuses on excellent performance: high resolution images, broadband and multi-frequency transducers, abundant obstetrics/gynecology software packages, dual USB ports, 256-frame cine loop, 56MB build-in image storage and various storage forms, etc., all of which make clinical diagnosis more functional and convenient.

DUS 8

Digital Ultrasonic Diagnostic Imaging System

Technical Specifications

General

Imaging mode: B,B+B,4B, B+M,M Gray scales: 256 Display: 14" non-interlaced Transducer frequency: 2.0 ~ 10MHz Transducer connector: 2 standard Beam-forming: Digital Beam-forming **Dynamic Receiving Focusing** Real-time Dynamic Aperture **Dynamic Frequency Scanning** Dynamic Apodization Tissue Harmonic Imaging **Tissue Specific Imaging** Scanning angle: from 30 to 155 degree (depending on transducers) Scanning depth (mm): from 20 to 250 (depending on transducers)

Imaging Processing:

Pre-processing: Dynamic range Edge enhancement Frame correlation Line correlation Smooth AGC 8-segment TGC adjustment IP (Image Process) Post-processing: Gray map Gamma correction

Rejection Left-right reverse Up-down reverse

Functions:

Cine loop: 256 frames bidirectional cine-loop Zoom: X1.0, X1.2, X1.4, X1.6, X2.0, X2.4, X3.0, X4.0 in distance Storage media: Built-in Flash, External USB-Memory stick Storage: 56MB permanent image Body mark: > 80 types Transducer auto-detection 16-segment acoustic power output adjustment

Measurement & Calculation:

B-mode: distance, circumference, area, volume, angle, Ratio,%stenosis M-mode: distance, time, velocity, heart rate (2 cycles), slope Software packages: abdomen, gynecology, obstetrics, urology, small parts, cardiology, orthopedics





All rights reserved.

Edan Instruments, Inc.

Date, Time, Probe Name, Probe Frequency, Frame Rate, Patient Name, Patient ID, Hospital Name, Measurement Values, Body Marks, Annotation, Probe Position, Full-image-region edit

ar array : L743

Endorectal: E743

Others:

Display

Endovaginal: E613

Peripheral port: Video output 1 VGA output port 1 USB port 2 DICOM3.0 1 (optional) 100V-240V ~ 50Hz/60Hz Power supply: 530mm(W) X 700mm(L) X1300mm(H) **Dimensions:** Net weight: 33ka

Multi-frequency transducers

Standard Configurations:

DUS 8 main unit 14" non-interlaced monitor Two transducer connectors 256 frames cine loop memory 56 MB built-in image storage Two USB ports Measurement & calculation software packages Convex array transducer: C363-1

Options:

Linear array transducer:	L743
Endorectal transducer:	E743
Endovaginal transducer:	E613
Micro-convex array transducer:	C321
Convex array transducer:	C343-1

Video printer

Laser printer Biopsy guide DICOM3.0 Footswitch





ENG-US-DUS8-V1.2-20091130